

SDCCAG8 Antibody



PACO12071

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SDCCAG8

Uniprot

Q86SQ7

Synonyms:

serologically defined colon cancer antigen 8; SDCCAG8; CCCAP; HSPC085; NY-CO-8 ;

Immunogen:

Human SDCCAG8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A