

# SDR39U1 Antibody



PACO12077

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SDR39U1

**Uniprot**

Q9NRG7

**Synonyms:**

short chain dehydrogenase/reductase family 39U, member 1; SDR39U1; C14orf124; HCDI ;

**Immunogen:**

Human SDR39U1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A