

SEC13 Antibody



PACO12080

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SEC13

Uniprot

P55735

Synonyms:

SEC13 homolog (*S. cerevisiae*); SEC13; D3S1231E; SEC13L1; SEC13R; npp-20 ;

Immunogen:

Human SEC13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A