

SEC63 Antibody



PACO12089

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SEC63

Uniprot

Q9UGP8

Synonyms:

SEC63 homolog (*S. cerevisiae*); SEC63; ERdj2; PRO2507; SEC63L ;

Immunogen:

Human SEC63.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A