

SERF2 Antibody



PACO12101

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SERF2

Uniprot

P84101

Synonyms:

small EDRK-rich factor 2; SERF2;4F5REL; FAM2C; FLJ20431; FLJ37527; FLJ38557; H4F5rel; HsT17089; MGC48826 ;

Immunogen:

Human SERF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A