

SERP1 Antibody



PACO12104

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SERP1

Uniprot

Q9Y6X1

Synonyms:

stress-associated endoplasmic reticulum protein 1; SERP1; MGC117327; MGC133321; MGC133322; RAMP4 ;

Immunogen:

Human SERP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A