

# SERPINB10 Antibody



PACO12110

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SERPINB10

**Uniprot**

P48595

**Synonyms:**

serpin peptidase inhibitor, clade B (ovalbumin), member 10; SERPINB10; PI10; bomapin ;

**Immunogen:**

Human SERPINB10.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A