

SRSF5 Antibody



PACO12149

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SRSF5

Uniprot

Q13243

Synonyms:

splicing factor, arginine/serine-rich 5; SFRS5; HRS; SRP40 ;

Immunogen:

Human SFRS5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A