

SGSM1 Antibody



PACO12170

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SGSM1

Uniprot

Q2NKQ1

Synonyms:

small G protein signaling modulator 1; SGSM1; KIAA1941; MGC133298; RUTBC2 ;

Immunogen:

Human SGSM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A