

SGSM3 Antibody



PACO12171

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SGSM3

Uniprot

Q96HU1

Synonyms:

small G protein signaling modulator 3; SGSM3; DKFZp761D051; MAP; RUSC3; RUTBC3 ;

Immunogen:

Human SGSM3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A