

SH2D4A Antibody



PACO12178

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SH2D4A

Uniprot

Q9H788

Synonyms:

SH2 domain containing 4A; SH2D4A; FLJ20967; SH2A ;

Immunogen:

Human SH2D4A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A