

SH3BP1 Antibody



PACO12181

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SH3BP1

Uniprot

Q9Y3L3

Synonyms:

SH3-domain binding protein 1;

Immunogen:

Human SH3BP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A