

SH3BP4 Antibody



PACO12182

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SH3BP4

Uniprot

Q9P0V3

Synonyms:

SH3-domain binding protein 4; SH3BP4; BOG25; TTP ;

Immunogen:

Human SH3BP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A