

# SH3BP5 Antibody



PACO12183

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SH3BP5

**Uniprot**

O60239

**Synonyms:**

SH3-domain binding protein 5 (BTK-associated); SH3BP5; SAB ;

**Immunogen:**

Human SH3BP5.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A