

# SH3D19 Antibody



PACO12184

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SH3D19

**Uniprot**

Q5HYK7

**Synonyms:**

SH3 domain containing 19; SH3D19; EBP; EVE1; Kryn; MGC105136; MGC118910; MGC118911; MGC118912; MGC118913; SH3P19 ;

**Immunogen:**

Human SH3D19.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A