

SHKBP1 Antibody



PACO12199

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SHKBP1

Uniprot

Q8TBC3

Synonyms:

SH3KBP1 binding protein 1; SHKBP1; PP203; Sb1 ;

Immunogen:

Human SHKBP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A