

SHROOM1 Antibody



PACO12204

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SHROOM1

Uniprot

Q2M3G4

Synonyms:

shroom family member 1; SHROOM1; APXL2; KIAA1960 ;

Immunogen:

Human SHROOM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A