

SLC12A4 Antibody



PACO12239

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC12A4

Uniprot

Q9UP95

Synonyms:

solute carrier family 12 (potassium/chloride transporters), member 4; SLC12A4; FLJ17069; FLJ40489; KCC1 ;

Immunogen:

Human SLC12A4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A