

SLC17A5 Antibody



PACO12244

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC17A5

Uniprot

Q9NRA2

Synonyms:

solute carrier family 17 (anion/sugar transporter), member 5; SLC17A5; AST; FLJ22227; FLJ23268; ISSD; NSD; SD; SIALIN; SIASD; SLD ;

Immunogen:

Human SLC17A5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A