

SLC22A18 Antibody



PACO12254

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC22A18

Uniprot

Q96BI1

Synonyms:

solute carrier family 22, member 18; SLC22A18; BWR1A; BWSCR1A; DKFZp667A184; HET; IMPT1; ITM; ORCTL2; SLC22A1L; TSSC5; p45-BWR1A ;

Immunogen:

Human SLC22A18.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A