

SLC25A2 Antibody



PACO12264

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SLC25A2

Uniprot

Q9BXI2

Synonyms:

solute carrier family 25 (mitochondrial carrier; ornithine transporter) member 2;
SLC25A2; MGC119151; MGC119153; ORC2; ORNT2 ;

Immunogen:

Human SLC25A2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A