

SLC25A25 Antibody



PACO12268

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC25A25

Uniprot

Q6KCM7

Synonyms:

solute carrier family 25 (mitochondrial carrier; phosphate carrier), member 25; SLC25A25; KIAA1896; MCSC; MGC105138; MGC119514; MGC119515; MGC119516; MGC119517; PCSCL; RP11-395P17.4; SCAMC-2 ;

Immunogen:

Human SLC25A25.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A