

# SLC25A3 Antibody



PACO12269

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC25A3

**Uniprot**

Q00325

**Synonyms:**

solute carrier family 25 (mitochondrial carrier; phosphate carrier), member 3; SLC25A3; OK/SW-cl.48; PHC ;

**Immunogen:**

Human SLC25A3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A