

SLC30A6 Antibody



PACO12281

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC30A6

Uniprot

Q6NXT4

Synonyms:

solute carrier family 30 (zinc transporter), member 6; SLC30A6; FLJ31101; MGC45055; MST103; MSTP103; ZNT6 ;

Immunogen:

Human SLC30A6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A