

# SLC30A7 Antibody



PACO12282

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SLC30A7

**Uniprot**

Q8NEW0

**Synonyms:**

solute carrier family 30 (zinc transporter), member 7; SLC30A7; DKFZp686M0368; ZNT7; ZnT-7; ZnTL2 ;

**Immunogen:**

Human SLC30A7.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A