

SLC38A3 Antibody



PACO12288

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SLC38A3

Uniprot

Q99624

Synonyms:

solute carrier family 38, member 3; SLC38A3; G17; NAT1; SN1 ;

Immunogen:

Human SLC38A3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A