

SLC39A8 Antibody



PACO12293

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SLC39A8

Uniprot

Q9C0K1

Synonyms:

solute carrier family 39 (zinc transporter), member 8; SLC39A8; BIGM103; LZT-Hs6; PP3105; ZIP8 ;

Immunogen:

Human SLC39A8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A