

# SLC43A1 Antibody



PACO12296

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

SLC43A1

**Uniprot**

O75387

**Synonyms:**

61kd;

**Immunogen:**

Human SLC43A1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.320&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A