

SLC45A2 Antibody



PACO12298

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC45A2

Uniprot

Q9UMX9

Synonyms:

solute carrier family 45, member 2; SLC45A2;1A1; AIM1; MATP; SHEP5 ;

Immunogen:

Human SLC45A2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A