

SLC4A8 Antibody



PACO12301

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC4A8

Uniprot

Q2Y0W8

Synonyms:

solute carrier family 4, sodium bicarbonate cotransporter, member 8; SLC4A8; DKFZp761B2318; FLJ46462; NBC3 ;

Immunogen:

Human SLC4A8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A