

SLC7A4 Antibody



PACO12306

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SLC7A4

Uniprot

O43246

Synonyms:

solute carrier family 7 (cationic amino acid, transporter, y+ system), member 4; SLC7A4; CAT-4; CAT4; HCAT3; MGC129976; MGC129977; VH ;

Immunogen:

Human SLC7A4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A