

SLC9C2 Antibody



PACO12308

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SLC9C2

Uniprot

Q5TAH2

Synonyms:

solute carrier family 9, member 11; SLC9A11; MGC43026; RP3-436N22.2 ;

Immunogen:

Human SLC9A11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A