

SMC6 Antibody



PACO12341

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SMC6

Uniprot

Q96SB8

Synonyms:

structural maintenance of chromosomes 6; SMC6; FLJ22116; FLJ35534; SMC6L1 ;

Immunogen:

Human SMC6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A