

SMTNL2 Antibody



PACO12349

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SMTNL2

Uniprot

Q2TAL5

Synonyms:

smoothelin-like 2; SMTNL2; FLJ42461; MGC131847; MGC138382 ;

Immunogen:

Human SMTNL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A