

SMU1 Antibody



PACO12350

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SMU1

Uniprot

Q2TAY7

Synonyms:

smu-1 suppressor of mec-8 and unc-52 homolog (C. elegans); SMU1; BWD; DKFZp761L0916; FLJ10805; FLJ10870; FLJ11970; MGC117363; RP11-54K16.3; SMU-1; fSAP57 ;

Immunogen:

Human SMU1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A