

# SP100 Antibody



PACO12405

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SP100

**Uniprot**

P23497

**Synonyms:**

SP100 nuclear antigen; SP100; DKFZp686E07254; FLJ00340; FLJ34579 ;

**Immunogen:**

Human SP100.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A