

SPACA3 Antibody



PACO12409

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SPACA3

Uniprot

Q8IXA5

Synonyms:

sperm acrosome associated 3; SPACA3;1700025M08Rik; ALLP17; CT54; LYC3; LYZL3; MGC119058; SLLP1 ;

Immunogen:

Human SPACA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A