

SPAG11A Antibody



PACO12410

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SPAG11A

Uniprot

Q6PDA7

Synonyms:

sperm associated antigen 11A; SPAG11A; HE2 ;

Immunogen:

Human SPAG11A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A