

# SPAG7 Antibody



PACO12415

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

SPAG7

**Uniprot**

O75391

**Synonyms:**

sperm associated antigen 7; SPAG7; ACRP; FSA-1; MGC20134 ;

**Immunogen:**

Human SPAG7.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A