

SPANXE Antibody



PACO12419

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SPANXE

Uniprot

Q8TAD1

Synonyms:

SPANX family, member E; SPANXE; MGC71908; SPANX-A2; SPANX-C; SPANX-D; SPANXA2; SPANXD ;

Immunogen:

Human SPANXE.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A