

SPATA13 Antibody



PACO12422

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SPATA13

Uniprot

Q96N96

Synonyms:

spermatogenesis associated 13; SPATA13; ASEF2; FLJ31208; FLJ35435; MGC129988; MGC129989 ;

Immunogen:

Human SPATA13.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A