

# SPATA2 Antibody



PACO12424

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SPATA2

**Uniprot**

Q9UM82

**Synonyms:**

spermatogenesis associated 2; SPATA2; FLJ13167; KIAA0757; PD1; tamo ;

**Immunogen:**

Human SPATA2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A