

SPATA22 Antibody



PACO12425

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SPATA22

Uniprot

Q8NHS9

Synonyms:

spermatogenesis associated 22; SPATA22; NYD-SP20 ;

Immunogen:

Human SPATA22.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A