

# SPCS1 Antibody



PACO12427

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SPCS1

**Uniprot**

Q9Y6A9

**Synonyms:**

signal peptidase complex subunit 1 homolog (S. cerevisiae); SPCS1; HSPC033; SPC1; SPC12; YJR010C-A ;

**Immunogen:**

Human SPCS1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A