

PACO12429

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SPERT

Uniprot

Q8NA61

Synonyms:

spermatid associated; SPERT; CBY2; FLJ35810; NURIT ;

Immunogen:

Human SPERT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A