

SPG20 Antibody



PACO12430

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

SPART

Uniprot

Q8N0X7

Synonyms:

spastic paraplegia 20 (Troyer syndrome); SPG20; KIAA0610; SPARTIN; TAHCCP1 ;

Immunogen:

Human SPG20.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A