

SPHK1 Antibody



PACO12431

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SPHK1

Uniprot

Q9NYA1

Synonyms:

sphingosine kinase 1; SPHK1; SPHK ;

Immunogen:

Human SPHK1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A