

PACO12441

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SPR

Uniprot

P35270

Synonyms:

sepiapterin reductase (7,8-dihydrobiopterin: NADP+ oxidoreductase); SPR; SDR38C1 ;

Immunogen:

Human SPR.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A