

SPZ1 Antibody



PACO12448

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SPZ1

Uniprot

Q9BXG8

Synonyms:

spermatogenic leucine zipper 1; SPZ1; FLJ25709; NYD-TSP1 ;

Immunogen:

Human SPZ1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A