

# SQRDL Antibody



PACO12450

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SQOR

**Uniprot**

Q9Y6N5

**Synonyms:**

sulfide quinone reductase-like (yeast); SQRDL; CGI-44 ;

**Immunogen:**

Human SQRDL

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A